

Relevo™



The super absorbent dressing for heavily exuding wounds

- Absorbs more than traditional absorbent dressings
- Absorbs under compression and pressure
- Binds the fluid inside the dressing without leaking and dripping
- Contributes to fewer changes, time saving and better patient comfort
- Absorbs vertically and can reduce the risk of maceration on healthy skin area
- Sterile, single-packed in peel-pouch
- CE marked (Class IIb, according to MDD 93/42/EEC)

Relevo™



	Relevo Dressing 10 x 10 cm	Relevo Dressing 10 x 10 cm	Relevo Dressing 10 x 22 cm	Relevo Dressing 20 x 22 cm
Product description	Code F5001 I/10	Code F5001 I/50	Code F5002I	Code F5004I
Absorption capacity (Tap water)	120 ml	120 ml	300 ml	600 ml
Absorption capacity (Saline solution 0,9 %)	50 ml	50 ml	120 ml	240 ml
Packing Pieces/box	10	50	50	50
Packing Inner boxes/ transport box	45	9	9	6
EAN Code	7350031911428	7350031911466	7350031910056	7350031910070
Reorder Code	1788549	1788581	1770023	1770031

- Use:** Single use only
- Can be cut:** No
- Indications:** Heavily exuding wounds
- Contains:** Cellulose and super polymers with a cover of polypropylene
- Sterile:** Yes
- Disposal:** Normal waste. Combustion forms carbon dioxide and water.
If contaminated follow local regulations
- Latex Free:** Yes
- TGA approved:** Number: I44120 Sponsor: Reliance Medical Pty Ltd